

CMS - 5100



Features

- +Suitable for adult, pediatric and neonatal patient.
- +2.8" High brightness TFT LCD display.
- +High brightness LED display of NIBP, SPO₂ and pulse rate.
- +Memory: Built-in Flash Memory.
- +Up to 2000 groups NIBP data storage and SPO₂, data up to 48 hours.
- +Built-in rechargeable Lithium ion polymer battery.

Specifications

NIBP

Method Oscillometry.

Mode Manual/Auto/Continuous.

Measurement Range 25-260mmHg. Resolution 1mmHg.

NIBP Accuracy ±5mmHg.

Alarm Systolic/Diastolic/Map.

SpO₂

Measurement Range 0%~100%

Resolution 1%.

Accuracy ±2%(70%~100%), 0%~69% unspecified.

Alarm range 0~100%.

Pulse Rate

Measurement Range 10-270bpm

Resolution 1bpm.

Accuracy ±2bpm.

