

Medical Oxygen Concentrator

RMS Oxygen Concentrator



Respiricare		Oxygen Concentrator
Model		A05
Flow rate		0.5-5 Lit
Oxygen Concentration		93% <u>±</u> 3%
Oxygen Output Pressure		0.04-0.06Mpa
Nebulization		No
Oxygen Sensor		Yes
Display		Oxygen Purity, Flowrate, Time function
Alaram		High-Low Pressure Alaram, Low Oxygen Concentrator Alaram,Power off Alaram.
Other Functions		Timer Function, Regular Maintenince alert
Noise Level		≤45dB
Power Supply		220 <u>±</u> 22V , 50Hz <u>±</u> 1Hz
Dimensions		310mm*312mm*510mm
Weight		18Kgs.

