


BIPAP

RMS BIPAP T25



Mode->S/T	
Mode	S/T
IPAP	4.0 cm-H ₂ O
EPAP	0.40 cm-H ₂ O
Rise Time	0.1 sec
Rate	10.0 bpm
Insp Time	1.5 sec
I-trigger	6



Full Face Mask

Model	BPAP25
Modes	CPAPS,T,S/T
IPAP Range	4-25 Cm H ₂ O
EPAP Range	4-20 Cm H ₂ O
CPAP Range	4-20 Cm H ₂ O
I-Trigger	Auto,1-6
E-Trigger	Auto,1-6
IPAP Rise Time	Auto,0.1-0.6s
Breath Rate	2-40 Bpm
Noise Level (db)	<28db
Ramp Time(minute)	0-60
Display Data	Ipap, Epap, Rate(bpm), Le, Vt, Vleak, Mvent,triggeri/e
Alarms	Tube Alarm, Apnea Alarm, High/low Rate Alarms & Power Failure
Dimensions	18'.27*11cm
Weight	1.7 Kg
Screen Size	4.3 Inch
Heated Humidifier	Integrated
Memory	Micro Sd Card - 8 Gb



We support service...