

# Syringe Pump

## HAWKMED - HK - 400



### Syringe Pump Features:

- History records 1600 record
- Unique intelligent self-recognition technology of syringes
- KVO, bolus and anti-bolus
- Applicable for, 10,20,30,50(60) ml syringes
- Compatible for various brands of syringes
- Various infusion modes to meet different clinical requirements
- Large colorful lcd screen
- 3 adjustable occlusion levels, dynamically display the pressure
- Night mode
- Double CPU ensuring safe infusion
- **Infusion mode:** Rate mode | Time mode | Body weight mode  
Intermittent infusion mode (optional)  
Drug library mode (optional)
- **Battery operating time:** About 7 hours at 25ml/h

