



## HAWKMED - HK - 100 II

### Infusion Pump Features

- History records
- Various infusion modes to meet different Clinical requirements
- KVO, Bolus and Anti-Bolus
- Drop sensor (optional)
- Real time display
- Free-flow protection
- Double CPU ensuring safe infusion
- **Infusion Rate:** 0.1-1200ml/h
- **Infusion Mode:** Rate Mode | Drip Mode | Time Mode | Body Weight Mode  
Drug Library Mode(optional) | Intermittent Infusion Mode(optional)
- **Purge Rate:** 600ml/h
- **Battery Operating Time:** About 6 Hours At 25ml/h

