

## BMC Auto CPAP



Modes	RESmart GII-Auto CPAP
Cover Color	Black
Lcd Screen	3.5 inches
Mode	CPAP/Auto
Pressure Range	4-20cm H <sub>2</sub> O
Pressure Increment	0.5cm H <sub>2</sub> O
Reslex	Off/1/2/3
Ramp	0-60 minutes
Sound	<30dBA
Alarms	Power Fail/High Leak
Humidity Level	Off /1-5
Auto On/off	Off(Deafault)
Data Storage	SD Card



**Full Face Mask**

