

BIPAP



BMC GII BPAP - Y-25T (With AVPAS)

Bi-level Positive Airway Pressure
With Ti Min/ Ti Max

3.5 inch Colorful LCD Display
Real time pressure & flow wave display.

Parameter Product	Y-25T
Mode	CPAP, S, T, S/T
Pressure	4-25cm H ₂ O
LCD Sreen	3.5inches
Weight	5.5 lb(2.5 kg) with humidifer
Dimensions(LxWxH)	11.4" x7.7"x5.3" with humidifer
Oximeter	optional
Data Storage	SD Card
Reslex	Off/1/2/3
Ramp P	4cm h ₂ o
Max Ramp	0-60 minutes
Ramp Time	0 - max ramp
Humidity Level	off/1-5
Ti Min	0.1s - ti max
Ti Max	0.3s - 3.0s



Full Face Mask

