

PULSE OXIMETER

Hand Held Pulse Oximeter

General Features:

- Bright, Easy-to-Read Large LED display
- Durable , Compact and Light Weight
- Audiable/ Visual Alaram with adjustable alaram limits
- Continoues Monitoring and Spot-Check
- Four Standard "AA" Size Rechanrgble batteries with charger Dock
- Suitable for Adult, Paediatric and Neonetal
- Nickel-Cadmium battery of 4.8V
- Battery backup: Approx. 12 hours
- Alaram indicator for SpO2, PR & Sensor disconnectingproblem with adjustable alaram limits
- Compatible with Nellcor SpO2 probe
- Auto Shut-off function when the finger is out



General Features:

- Measurement range for SpO2 : 0-100%
- Accuracy for SpO2/PR : +2%
- Measurement range for PR: 30-250 bpm
- Alaram silence interval range: 15,30,60,90,120sec.

Finger Tip Pulse Oximeter

• Model: FTPCMS50D

• General Features:

- Colour Display with High Brightness
- Standard 2 No's .AAA battery size Alkaline Battery *1.5V each
- Audio and visual alaram for Spo2 and Pulse rate with Beep
- Auto shut-off function when the finger is out
- 4 Direction display with six level brightness adjustment
- Battery backup – Approx.25 hours

